



### **Product description**

Order number	97.8100.111
Product description	Pressure infusion cuff Standard, up to 500 ml, fixed insertion, latex-free

### **Product characteristics**

Type of cap	fixed insertion
Application	Transfusion and infusion

### **Size**

Max. volume	500 ml
Width of product	16 cm
Height of product	8 cm
Length of product	28 cm

### **Material & colours**

Colour of product	red
-------------------	-----

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product in accordance with MDD or IVDD	Medical device
CE certified	CE
Latex-free	yes
Lot or serial number	no

### **Packaging**

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0039 cbm
Weight of product	0,25 kg
Weight of case	0,42 kg
EAN of inner box	4028339311020
EAN case	4028339311013

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

