



Product description

Order number	97.8100.210
Product description	Cover, Replacement part for FR pressure infusion cuff

Product characteristics

Application	Transfusion and infusion
Suitable for	Pressure infusion cuff

Size

Max. volume	1000 ml
Width of product	14.5 cm
Height of product	2 cm
Length of product	27 cm

Material & colours

Colour of product	red
-------------------	-----

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 5
Case volume 0.0039 cbm
Weight of product 0.16 kg
Weight of case 0.97 kg
EAN of inner box 4028339311228
EAN case 4028339311211

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

