



### **Product description**

Order number	97.8100.410
Product description	Disposable pressure infusion cuff up to 500 ml, fixed insertion

### **Product characteristics**

Type of cap	fixed insertion
Application	Transfusion and infusion

### **Size**

Max. volume	500 ml
Width of product	15 cm
Height of product	6 cm
Length of product	24 cm

### **Material & colours**

Colour of product	grey
-------------------	------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product in accordance with MDD or IVDD	Medical device
CE certified	CE
Lot or serial number	yes

### **Packaging**

Minimum order qty.	2
Piece(s) / inner box	2
Piece(s) / case	20
Piece(s) / pallet	1440
Case volume	0,0183 cbm
Weight of product	0,1295 kg
Weight of case	3,6 kg
EAN of inner box	4028339311426
EAN case	4028339311419

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

