



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 97.8100.410   |
| Product description | Disposable pressure infusion cuff up to 500 ml, fixed insertion, with manometer |

### **Product characteristics**

|             |                          |
|-------------|--------------------------|
| Type of cap | fixed insertion          |
| Application | Transfusion and infusion |

### **Size**

|                   |        |
|-------------------|--------|
| Max. volume       | 500 ml |
| Width of product  | 15 cm  |
| Height of product | 6 cm   |
| Length of product | 24 cm  |

### **Material & colours**

|                   |      |
|-------------------|------|
| Colour of product | grey |
|-------------------|------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

|                  |                |
|------------------|----------------|
| Product category | Medical device |
| CE certified     | CE             |
| Batched          | yes            |

### **Packaging**

|                      |               |
|----------------------|---------------|
| Minimum order qty.   | 10            |
| Piece(s) / inner box | 10            |
| Piece(s) / case      | 20            |
| Piece(s) / pallet    | 1440          |
| Case volume          | 0,0183 cbm    |
| Weight of product    | 0,1295 kg     |
| Weight of case       | 3,6 kg        |
| EAN of inner box     | 4028339311426 |
| EAN case             | 4028339311419 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

