



Product description

Order number	79.8200.000
Product description	TRANSMED 10 Transfusion device, with integrated cascade filter, 175 µm, 40 µm, 10 µm, spike colour: yellow

Product characteristics

Application	Transfusion and infusion
-------------	--------------------------

Size

Length of tube	150 cm
----------------	--------

Material & colours

Colour of product	yellow
-------------------	--------

Purity & certification

Product category	Medical device
CE certified	CE 0482 (MedCert)
Purity standard	sterile
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	Trioctyl trimellitate (TOTM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Batched yes

Packaging

Minimum order qty.	100
Type of inner packaging	blister
Piece(s) / inner box	100
Piece(s) / case	100
Case volume	0.0566 cbm
Weight of product	0.0492 kg
Weight of case	5.95 kg
EAN of inner box	4028339420128
EAN case	4028339420111

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

