



Product description

Order number	79.8200.100
Product description	Transfusion device, with integrated cascade filter, 175 µm, 40 µm, spike colour: red

Product characteristics

Application	Transfusion and infusion
Connection	Spike

Size

Length of tube	150 cm
----------------	--------

Material & colours

Colour of product	red
-------------------	-----

Purity & certification

Product category	Medical device
CE certified	CE0482 (Medcert)
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	Trioctyl trimellitate (TOTM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Batched yes

Packaging

Minimum order qty.	100
Piece(s) / inner box	100
Piece(s) / case	100
Case volume	0,0566 cbm
Weight of product	0,05 kg
Weight of case	6 kg
EAN of inner box	4028339420524
EAN case	4028339420517

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

