



### **Product description**

Order number	79.8200.440
Product description	Transfusion device, with integrated filter, 40 µm, spike colour: green

### **Product characteristics**

Application	Transfusion and infusion
Filter	yes
Pore size	40

### **Size**

Length of tube	150 cm
----------------	--------

### **Material & colours**

Colour of product	green
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category	Medical device
CE certified	CE0482 (Medcert)
Purity standard	sterile
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	Trioctyl trimellitate (TOTM)
Batched	yes

### **Packaging**

Minimum order qty.	100
Piece(s) / inner box	100
Piece(s) / case	100
Piece(s) / pallet	2800
Case volume	0,0509 cbm
Weight of product	0,0405 kg
Weight of case	5,05 kg
EAN of inner box	4028339420722
EAN case	4028339420715

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

