



Product description

Order number	79.8089.219
Product description	CSE set, set components: REGANESTH® pencil-point spinal needle 27G x 125 mm, Tuohy needle 18G x 90 mm, epidural catheter 21G x 900 mm, Tuohy-Borst adapter, green hub insert, flat filter 0.2 µm, L.O.R. LL syringe

Product characteristics

Application	regional anaesthesia
Filter	yes
Pore size	0,2
Point sharpening	Tuohy, pencil-point
Connection	Luer

Size

Diameter of needle (external)	18 G
Length of needle	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	Medical device
CE certified	CE supplier > no Sarstedt CE
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	5
Piece(s) / inner box	5
Piece(s) / case	50
Depth of case	464 mm
Width of case	313 mm
Height of case	220 mm
Case volume	0,032 cbm
Weight of product	0,0226 kg
Weight of case	2,34 kg
EAN of inner box	4028339180824
EAN case	4028339180817

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

