



## Product description

Order number	79.8096.185
Product description	Plexus catheter set, set components: REGANESTH® unoplex 18G x 50 mm, plexus catheter 21G x 400 mm with mandrin, injection tube 60 cm, Tuohy-Borst adapter, hub insert, flat filter 0.2 µm with fixation

## Product characteristics

Application	regional anaesthesia
Filter	yes
Pore size	0,2
Point sharpening	Facet
Connection	Luer

## Size

Diameter of needle (external)	18 G / 1,35 mm
Length of needle	50 mm
Length of tube	60 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product in accordance with MDD or IVDD	Medical device
CE certified	CE supplier > no Sarstedt CE
Sterilisation	Gas sterilisation
Plasticiser	Trioctyl trimellitate (TOTM)
Lot or serial number	yes

### **Packaging**

Minimum order qty.	5
Piece(s) / inner box	5
Piece(s) / case	80
Piece(s) / pallet	3840
Depth of box	203 mm
Width of box	151 mm
Height of box	56 mm
Depth of case	464 mm
Width of case	313 mm
Height of case	220 mm
Case volume	0,032 cbm
Weight of product	0,04 kg
Weight of case	4,34 kg
EAN of inner box	4028339220421
EAN case	4028339220414

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

