



### Product description

Order number	79.8096.285
Product description	Plexus catheter set, set components: REGANESTH® unoplex 18G x 50 mm, plexus catheter 21G x 500 mm, central opening, with mandrin and Luer hub insert, injection tube 60 cm, Tuohy-Borst adapter, flat filter 0.2 µm with fixation

### Product characteristics

Application	regional anaesthesia
Filter	yes
Pore size	0,2
Point sharpening	Facet
Connection	Luer

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Diameter of needle (external)	18 G
Length of needle	50 mm
Length of tube	60 cm

### **Purity & certification**

Product in accordance with MDD or IVDD	Medical device
CE certified	CE supplier > no Sarstedt CE
Sterilisation	Gas sterilisation
Plasticiser	Trioctyl trimellitate (TOTM)
Lot or serial number	yes

### **Packaging**

Minimum order qty.	5
Piece(s) / inner box	5
Piece(s) / case	80
Depth of case	464 mm
Width of case	313 mm
Height of case	220 mm
Case volume	0,032 cbm
Weight of product	0,0261 kg
Weight of case	3,69 kg
EAN of inner box	4028339220520
EAN case	4028339220513

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

