



### **Product description**

Order number	79.8096.050
Product description	Plexus needle, 22G x 50 mm, can be visualised with ultrasound, connection: Luer (f), injection tube detachable, tube length: 60 cm

### **Product characteristics**

Application	regional anaesthesia
Point sharpening	Facet
Connection	Luer (f)

### **Size**

Diameter of needle (external)	22 G / 0,7 mm
Length of needle	50 mm
Length of tube	60 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product in accordance with MDD or IVDD	Medical device
CE certified	CE0482 (Medcert)
Sterilisation	Gas sterilisation
Plasticiser	Trioctyl trimellitate (TOTM)
Lot or serial number	yes

### **Packaging**

Minimum order qty.	10
Piece(s) / inner box	10
Piece(s) / case	160
Piece(s) / pallet	5760
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of product	0,0107 kg
Weight of case	2,93 kg
EAN of inner box	4028339210620
EAN case	4028339210613

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

