



Product description

Order number	79.8096.180
Product description	Plexus needle, 18G x 80 mm, can be visualised with ultrasound, connection: Luer

Product characteristics

Application	regional anaesthesia
Point sharpening	Facet
Connection	Luer

Size

Diameter of needle (external)	18 G / 1,35 mm
Length of needle	80 mm

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE0482 (Medcert)
Sterilisation	Gas sterilisation
Lot or serial number	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	20
Piece(s) / inner box	20
Piece(s) / case	400
Piece(s) / pallet	19200
Depth of box	203 mm
Width of box	111 mm
Height of box	56 mm
Depth of case	464 mm
Width of case	313 mm
Height of case	220 mm
Case volume	0,032 cbm
Weight of product	0,0094 kg
Weight of case	4,84 kg
EAN of inner box	4028339211122
EAN case	4028339211115

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

