



Product description

Order number	79.8096.044
Product description	Plexus needle, 25G x 40 mm, can be visualised with ultrasound, connection: NRFit/ISO 80369-6, injection tube detachable, tube length: 60 cm

Product characteristics

Application	regional anaesthesia
Point sharpening	Facet
Connection	NRFit/ISO 80369-6

Size

Diameter of needle (external)	25 G
Length of needle	40 mm
Length of tube	60 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	Medical device
CE certified	CE0482 (Medcert)
Pyrogen-free	no

Packaging

Piece(s) / inner box	10
Piece(s) / case	160
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0 cbm
Weight of case	1,75 kg
EAN of inner box	4028339211528
EAN case	4028339211511

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

