



### **Product description**

Order number	79.8096.054
Product description	REGANESTH® Plexus needle unoplex NRFit, 22G x 50 mm, can be visualised with ultrasound, connection: NRFit/ISO 80369-6, injection tube detachable, tube length: 60 cm

### **Product characteristics**

Application	regional anaesthesia
Point sharpening	Facet
Connection	NRFit/ISO 80369-6

### **Size**

Diameter of needle (external)	22 G
Length of needle	50 mm
Length of tube	60 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**REGANESTH® Plexus needle unoplex NRFit 22G**  
**x 50 mm**

**Page 2**

**Purity & certification**

Product category	Medical device
CE certified	CE 0482 (MedCert)
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	10
Piece(s) / inner box	10
Piece(s) / case	160
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of product	0,0047 kg
Weight of case	2,58 kg
EAN of inner box	4028339210927
EAN case	4028339210910

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-01-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**