



Product description

| | |
|---------------------|--|
| Order number | 79.8096.104 |
| Product description | Plexus needle, 20G x 100 mm, can be visualised with ultrasound, connection: NRFit/ISO 80369-6, injection tube detachable, tube length: 60 cm |

Product characteristics

| | |
|------------------|----------------------|
| Application | regional anaesthesia |
| Point sharpening | Facet |
| Connection | NRFit/ISO 80369-6 |

Size

| | |
|-------------------------------|--------|
| Diameter of needle (external) | 20 G |
| Length of needle | 100 mm |
| Length of tube | 60 cm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
REGANESTH® Plexus needle unoplex NRFit 20G
x 100 mm

Page 2

Purity & certification

| | |
|------------------|------------------|
| Product category | Medical device |
| CE certified | CE0482 (Medcert) |
| Pyrogen-free | no |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 10 |
| Piece(s) / case | 160 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 300 mm |
| Case volume | 0 cbm |
| Weight of case | 1,75 kg |
| EAN of inner box | 4028339211627 |
| EAN case | 4028339211610 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT