



### **Product description**

Order number	79.8086.000
Product description	Flat filter, pore size: 0.2 µm, bacterial filter for in-line filtration during the injection/infusion of local anaesthetics

### **Product characteristics**

Application	regional anaesthesia
Filter	yes
Pore size	0,2
Connection	Luer

### **Size**

Width of product	30 mm
Height of product	7,2 mm
Length of product	62,8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product in accordance with MDD or IVDD	Medical device
CE certified	CE0482 (Medcert)
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Lot or serial number	yes

### **Packaging**

Minimum order qty.	20
Piece(s) / inner box	20
Piece(s) / case	320
Piece(s) / pallet	15360
Depth of box	203 mm
Width of box	151 mm
Height of box	56 mm
Depth of case	464 mm
Width of case	313 mm
Height of case	220 mm
Case volume	0,032 cbm
Weight of product	0,01 kg
Weight of case	4,33 kg
EAN of inner box	4028339190328
EAN case	4028339190311

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

