



Product description

Order number	79.8081.803
Product description	REGANESTH® Spinal needle, pencil-point, 25G x 103 mm, connection: Luer (f), with introducer: 20G x 30 mm

Product characteristics

Application	regional anaesthesia
Coding colour	orange
Introducer	20 G / 0,9 mm
Point sharpening	pencil-point
Connection	Luer (f)

Size

Diameter of needle (external)	25 G / 0,5 mm
Length of needle	103 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	Medical device
CE certified	CE 0482 (MedCert)
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	20
Type of inner packaging	blister
Piece(s) / inner box	20
Piece(s) / case	400
Depth of box	203 mm
Width of box	111 mm
Height of box	56 mm
Depth of case	464 mm
Width of case	313 mm
Height of case	220 mm
Case volume	0.032 cbm
Weight of product	0.0052 kg
Weight of case	3.18 kg
EAN of inner box	4028339111828
EAN case	4028339111811

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

