



### **Product description**

Order number	79.8081.834
Product description	Spinal needle, pencil-point, 24G x 150 mm, connection: NRFit/ISO 80369-6, with introducer: 20G x 40 mm

### **Product characteristics**

Application	regional anaesthesia
Coding colour	violet
Introducer	20 G / 0,9 mm
Point sharpening	pencil-point
Connection	NRFit/ISO 80369-6

### **Size**

Diameter of needle (external)	24 G / 0,55 mm
Length of needle	150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**REGANESTH® Spinal needle pencil-point NRFit**  
**24G x 150 mm with introducer**

**Page 2**

**Purity & certification**

Product category	Medical device
CE certified	CE0482 (Medcert)
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	20
Piece(s) / inner box	20
Piece(s) / case	320
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of product	0,0035 kg
Weight of case	2,96 kg
EAN of inner box	4028339113921
EAN case	4028339113914

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**