



Product description

Order number 79.8082.024
Product description Spinal needle, Quincke, 22G x 38 mm, connection:
NRFit/ISO 80369-6

Product characteristics

Application regional anaesthesia
Coding colour black
Point sharpening Quincke
Connection NRFit/ISO 80369-6

Size

Diameter of needle (external) 22 G / 0,7 mm
Length of needle 38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
REGANESTH® Spinal needle Quincke NRFit 22G
x 38 mm

Page 2

Purity & certification

Product category	Medical device
CE certified	CE0482 (Medcert)
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	20
Piece(s) / inner box	20
Piece(s) / case	400
Depth of case	464 mm
Width of case	313 mm
Height of case	220 mm
Case volume	0,032 cbm
Weight of product	0,0033 kg
Weight of case	3,2 kg
EAN of inner box	4028339121421
EAN case	4028339121414

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT