



### **Product description**

Order number	79.8082.184
Product description	REGANESTH® Spinal needle, Quincke, 22G x 180 mm, connection: NRFit/ISO 80369-6

### **Product characteristics**

Application	regional anaesthesia
Coding colour	black
Point sharpening	Quincke
Connection	NRFit/ISO 80369-6

### **Size**

Diameter of needle (external)	22 G / 0,7 mm
Length of needle	180 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**REGANESTH® Spinal needle Quincke NRFit 22G**  
**x 180 mm**

**Page 2**

**Purity & certification**

Product category	Medical device
CE certified	CE 0482 (MedCert)
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	20
Type of inner packaging	blister
Piece(s) / inner box	20
Piece(s) / case	320
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0068 kg
Weight of case	4.02 kg
EAN of inner box	4028339121629
EAN case	4028339121612

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-07-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**