



Product description

Order number	97.8300.501
Product description	Label dispenser, size (LxWxH): 18 x 21 x 14 cm, tabletop or wall-mounted, holds 10 rolls of syringe labels or 4 rolls of perfusor labels

Product characteristics

Application	Labelling
-------------	-----------

Size

Width of product	21 cm
Height of product	14 cm
Length of product	18 cm

Material & colours

Colour of product	white
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 165
Case volume 0.0088 cbm
Weight of product 2.06 kg
Weight of case 2.23 kg
EAN of inner box 4028339550221
EAN case 4028339550214

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

