



Product description

Order number	74.5766.001
Product description	Velcro band, size: S (380?430 mm), material: nylon
Aid number	15.25.11.0039
PZN number	8406664

Size

Length of product	430 mm
-------------------	--------

Material & colours

Product material	nylon
Colour of product	white

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE
Lot or serial number	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	10
Piece(s) / inner box	1
Piece(s) / case	10
Piece(s) / pallet	3240
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0,0051 cbm
Weight of product	0,016 kg
Weight of case	0,29 kg
EAN of inner box	4038917136600
EAN case	4038917027205

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

