



Product description

Order number	74.5767.001
Product description	UroStretch, size: S, material: PE
Aid number	15.25.11.1034
PZN number	8406948

Material & colours

Product material	polyester (PE)
Colour of product	white

Purity & certification

Product category	Medical device
CE certified	CE supplier
Batched	yes

Packaging

Minimum order qty.	10
Type of inner packaging	case
Piece(s) / inner box	10
Piece(s) / case	10
Piece(s) / pallet	2340

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
UroStretch, size: S

Page 2

Depth of case	283 mm
Width of case	173 mm
Height of case	135 mm
Case volume	0.0066 cbm
Weight of product	0.025 kg
Weight of case	0.37 kg
EAN of inner box	5703936100107
EAN case	5703936100145

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT