



Product description

| | |
|---------------------|--|
| Order number | 74.5767.002 |
| Product description | UroStretch, size: M, material: polyester |
| Aid number | 15.25.11.1034 |
| PZN number | 8406954 |

Material & colours

| | |
|-------------------|-----------|
| Product material | polyester |
| Colour of product | white |

Purity & certification

| | |
|--|------------------------------|
| Product in accordance with MDD or IVDD | Medical device |
| CE certified | CE supplier > no Sarstedt CE |
| Lot or serial number | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

| | |
|----------------------|---------------|
| Minimum order qty. | 10 |
| Piece(s) / inner box | 10 |
| Piece(s) / case | 10 |
| Piece(s) / pallet | 2340 |
| Depth of case | 283 mm |
| Width of case | 173 mm |
| Height of case | 135 mm |
| Case volume | 0,0066 cbm |
| Weight of product | 0,043 kg |
| Weight of case | 0,55 kg |
| EAN of inner box | 4038917136648 |
| EAN case | 4038917027243 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

