



Product description

Order number	74.5767.003
Product description	UroStretch, size: L, material: polyester
Aid number	15.25.11.1034
PZN number	8406960

Material & colours

Product material	polyester
Colour of product	white

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE supplier > no Sarstedt CE
Lot or serial number	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	10
Piece(s) / inner box	10
Piece(s) / case	10
Piece(s) / pallet	2340
Depth of case	283 mm
Width of case	173 mm
Height of case	135 mm
Case volume	0,0066 cbm
Weight of product	0,043 kg
Weight of case	0,55 kg
EAN of inner box	4038917136655
EAN case	4038917027250

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

