



### Product description

Order number	74.5218.001
Product description	Urine bag, 1.5 l, with anti-reflux valve, drainage valve, ETO treated, tube length: 90 cm
PZN number	16333672

### Product characteristics

Colour of print/label	blue
Design	drainage valve
Application	Urine and secretion drainage
Volume of work	1,5 l
Anti-reflux valve	yes

### Size

Sample volume	1,5 l
Max. volume	1,5 l
Diameter of tube	6,6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

### **Purity & certification**

Product in accordance with MDD or IVDD	Medical device
CE certified	CE
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Lot or serial number	yes

### **Packaging**

Minimum order qty.	150
Piece(s) / inner box	150
Piece(s) / case	150
Piece(s) / pallet	9600
Case volume	0,0279 cbm
Weight of product	0,0353 kg
Weight of case	5,8 kg
EAN of inner box	4038917388665
EAN case	4038917388658

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

