



Product description

Order number	74.5218.101
Product description	Urine bag, 1.5 l, with anti-reflux valve, drainage valve, tube length: 90 cm
PZN number	16333689

Product characteristics

Colour of print/label	blue
Design	drainage valve
Application	Urine and secretion drainage
Volume of work	1,5 l
Anti-reflux valve	yes

Size

Sample volume	1,5 l
Max. volume	1,5 l
Diameter of tube	6,6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Urine bag, 1.5 l, with anti-reflux valve, drainage valve

Page 2

Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE
Lot or serial number	yes

Packaging

Minimum order qty.	150
Piece(s) / inner box	150
Piece(s) / case	150
Piece(s) / pallet	9600
Depth of case	387 mm
Width of case	292 mm
Height of case	215 mm
Case volume	0,0243 cbm
Weight of product	0,0352 kg
Weight of case	5,8 kg
EAN of inner box	4038917354530
EAN case	4038917354523

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-07-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT