



Product description

Order number	74.5281.005
Product description	Secretion bag, 1.5 l, sterile, tube length: 90 cm, 1 piece(s)/blister
PZN number	16333614

Product characteristics

Colour of print/label	blue
Application	Drainage of exudate
Volume of work	1,5 l
Anti-reflux valve	no

Size

Sample volume	1,5 l
Max. volume	1,5 l
Diameter of tube	6,6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE 0124 (DEKRA)
Purity standard sterile
Sterilisation Gas sterilisation
Batched yes

Packaging

Minimum order qty. 100
Type of inner packaging blister
Piece(s) / inner box 100
Piece(s) / case 100
Piece(s) / pallet 4800
Depth of case 387 mm
Width of case 292 mm
Height of case 265 mm
Case volume 0,0299 cbm
Weight of product 0,036 kg
Weight of case 4,1 kg
EAN of inner box 4038917354561
EAN case 4038917354554

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

