



Product description

Order number	74.5281.101
Product description	Secretion bag, 1.5 l, tube length: 90 cm
PZN number	16333620

Product characteristics

Colour of print/label	blue
Application	Drainage of exudate
Volume of work	1,5 l
Anti-reflux valve	no

Size

Sample volume	1,5 l
Max. volume	1,5 l
Diameter of tube	6,6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE
Lot or serial number	yes

Packaging

Minimum order qty.	200
Piece(s) / inner box	200
Piece(s) / case	200
Piece(s) / pallet	9600
Depth of case	387 mm
Width of case	292 mm
Height of case	265 mm
Case volume	0,0299 cbm
Weight of product	0,0325 kg
Weight of case	7 kg
EAN of inner box	4038917354592
EAN case	4038917354585

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

