



Product description

Order number 74.5282.001
Product description Secretion bag, 2 l, ETO treated, tube length: 90 cm
PZN number 16333637

Product characteristics

Colour of print/label blue
Application Drainage of exudate
Volume of work 2 l
Anti-reflux valve no

Size

Sample volume 2 l
Max. volume 2 l
Diameter of tube 6,6 mm
Length of tube 90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Lot or serial number	yes

Packaging

Minimum order qty.	200
Piece(s) / inner box	200
Piece(s) / case	200
Piece(s) / pallet	9600
Case volume	0,0279 cbm
Weight of product	0,035 kg
Weight of case	7,4 kg
EAN of inner box	4038917388818
EAN case	4038917388801

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

