



Product description

Order number	74.5282.005
Product description	Secretion bag, 2 l, sterile, tube length: 90 cm, 1 piece(s)/blister
PZN number	16333643

Product characteristics

Colour of print/label	blue
Application	Drainage of exudate
Volume of work	2 l
Anti-reflux valve	no

Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	6,6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE0124 (Dekra)
Purity standard	sterile
Sterilisation	Gas sterilisation
Lot or serial number	yes

Packaging

Minimum order qty.	100
Type of inner packaging	blister
Piece(s) / inner box	1
Piece(s) / case	100
Piece(s) / pallet	9600
Case volume	0,0223 cbm
Weight of product	0,041 kg
Weight of case	4,4 kg
EAN of inner box	4038917388870
EAN case	4038917388863

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

