



Product description

Order number	74.5283.005
Product description	Urine bag, 1.5 l, with anti-reflux valve, sterile, tube length: 90 cm, 1 piece(s)/blister
PZN number	16333761

Product characteristics

Colour of print/label	blue
Application	Urine and secretion drainage
Volume of work	1,5 l
Anti-reflux valve	yes

Size

Sample volume	1,5 l
Max. volume	1,5 l
Diameter of tube	6,6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE0124 (Dekra)
Purity standard	sterile
Sterilisation	Gas sterilisation
Lot or serial number	yes

Packaging

Minimum order qty.	100
Type of inner packaging	blister
Piece(s) / inner box	100
Piece(s) / case	100
Piece(s) / pallet	4800
Depth of case	387 mm
Width of case	292 mm
Height of case	265 mm
Case volume	0,0299 cbm
Weight of product	0,0358 kg
Weight of case	4,1 kg
EAN of inner box	4038917354622
EAN case	4038917354615

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

