



Product description

Order number	74.5284.005
Product description	Urine bag, 2 l, with anti-reflux valve, sterile, tube length: 90 cm, 1 piece(s)/blister
PZN number	16335636

Product characteristics

Colour of print/label	blue
Application	Urine and secretion drainage
Volume of work	2 l
Anti-reflux valve	yes

Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE supplier
Purity standard sterile
Sterilisation Gas sterilisation
Batched yes

Packaging

Minimum order qty. 100
Type of inner packaging blister
Piece(s) / inner box 1
Piece(s) / case 100
Piece(s) / pallet 4800
Case volume 0.0223 cbm
Weight of product 0.041 kg
Weight of case 4.4 kg
EAN of inner box 6971243661034
EAN case 6971243660792

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

