



### **Product description**

Order number	74.5284.101
Product description	Urine bag, 2 l, with anti-reflux valve, tube length: 90 cm, 200 piece(s)/case
PZN number	16333749

### **Product characteristics**

Colour of print/label	blue
Application	Urine and secretion drainage
Volume of work	2 l
Anti-reflux valve	yes

### **Size**

Sample volume	2 l
Max. volume	2 l
Diameter of tube	6,6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

### **Purity & certification**

Product category	Medical device
CE certified	CE supplier
Batched	yes

### **Packaging**

Minimum order qty.	200
Type of inner packaging	case
Piece(s) / inner box	200
Piece(s) / case	200
Piece(s) / pallet	9600
Case volume	0,0279 cbm
Weight of product	0,035 kg
Weight of case	7,4 kg
EAN of inner box	4038917388993
EAN case	4038917388986

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

