




Product Specification
Gall secretion bag, 1.5 l, tube length: 90 cm,
sterile

Page 1



Product description

Order number	74.5286.005
Product description	Gall secretion bag, 1.5 l, sterile, tube length: 90 cm, 1 piece(s)/blister

Product characteristics

Colour of print/label	blue
Application	Drainage of exudate
Volume of work	1,5 l
Anti-reflux valve	no

Size

Sample volume	1,5 l
Max. volume	1,5 l
Diameter of tube	6,6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-07-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Gall secretion bag, 1.5 l, tube length: 90 cm,
sterile

Page 2

Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

Purity & certification

Product in accordance with MDD or IVDD Medical device
CE certified CE0124 (Dekra)
Purity standard sterile
Sterilisation Gas sterilisation
Lot or serial number yes

Packaging

Minimum order qty. 100
Type of inner packaging blister
Piece(s) / inner box 1
Piece(s) / case 100
Piece(s) / pallet 4800
Case volume 0,0299 cbm
Weight of product 0,038 kg
Weight of case 4,1 kg
EAN of inner box 4038917532709
EAN case 4038917532693

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-07-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT