



### **Product description**

Order number 97.8300.400  
Product description Fixing clamp, for fixation of the catheter quiver, height: 4 cm

### **Product characteristics**

Application Ward requirement

### **Size**

Width of product 4 cm  
Height of product 4 cm  
Length of product 9,5 cm

### **Material & colours**

Product material steel  
Colour of product black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category                   no medical device | IVD  
Batched                               no

### **Packaging**

Minimum order qty.               1  
Piece(s) / inner box               1  
Piece(s) / case                    1  
Case volume                       0,0039 cbm  
Weight of product                0,18 kg  
Weight of case                    0,35 kg  
EAN of inner box                 4028339560220  
EAN case                          4028339560213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

