



Product description

Order number	97.8300.450
Product description	Catheter quiver, for the storage of suction catheters, central-venous catheters etc., attachment: rail mounting bracket, height: 45 cm, Ø: 10 cm, removable base

Product characteristics

Application	Ward requirement
-------------	------------------

Size

Diameter	10 cm
Height of product	45 cm

Material & colours

Product material	Polypropylene (PP)
Colour of product	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 16
Case volume 0,02 cbm
Weight of product 0,0589 kg
Weight of case 0,95 kg
EAN of inner box 4028339560329
EAN case 4028339560312

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

