



Product description

Order number	97.8300.455
Product description	TRANSKAT Quiver base, Replacement part for TRANSKAT catheter quiver, Ø: 10 cm

Product characteristics

Application	Ward requirement
-------------	------------------

Size

Diameter	10 cm
----------	-------

Material & colours

Product material	Polypropylene (PP)
Colour of product	red

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0.0039 cbm
Weight of product	0.02 kg
Weight of case	0.19 kg
EAN of inner box	4028339560428
EAN case	4028339560411

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

