



Product description

Order number	97.8300.460
Product description	Protective cover, Cap for sealing the TRANSKAT catheter quiver, height: 25 cm, Ø: 10 cm

Product characteristics

Application	Ward requirement
-------------	------------------

Size

Diameter	10 cm
Height of product	25 cm

Material & colours

Product material	Polycarbonate (PC)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Case volume 0,012 cbm
Weight of product 0,23 kg
Weight of case 0,73 kg
EAN of inner box 4028339560527
EAN case 4028339560510

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

