



Product description

Order number	97.8300.470
Product description	Extension unit, for extending TRANSKAT catheter quiver by up to 11 cm, height: 14 cm, Ø: 10 cm

Product characteristics

Application	Ward requirement
-------------	------------------

Size

Diameter	10 cm
Height of product	14 cm

Material & colours

Colour of product	red
-------------------	-----

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0053 cbm
Weight of product	0,12 kg
Weight of case	0,42 kg
EAN of inner box	4028339561326
EAN case	4028339561319

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

