



### **Product description**

Order number	71.9924.515
Product description	Forceps, height: 11 mm, green, for Collection of diagnostic specimens

### **Product characteristics**

Application	Collection of diagnostic specimens
-------------	------------------------------------

### **Size**

Width of product	28 mm
Height of product	11 mm
Length of product	108 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Latex-free	yes
Pyrogen-free	no
Safety data sheet	nein

### **Packaging**

Minimum order qty.	3000
Piece(s) / inner box	3000
Piece(s) / case	3000
Piece(s) / pallet	90000
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0,0834 cbm
Weight of case	12 kg
EAN of inner box	4038917155281
EAN case	4038917398398

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

