



**Product Specification**  
**SURGISPONGE, impregnated with 4%  
chlorhexidine gluconate solution**

**Page 1**



**Product description**

Order number	54.9912.200
Product description	SURGISPONGE, Pre-operative scrubsponge for use as a total body wash for patients prior to surgery, impregnated with 4% chlorhexidine gluconate solution, material: PUR, (LxWxH): 210 x 63 x 45 mm

**Product characteristics**

Preparation concentration	4% chlorhexidine gluconate solution
Application	Pre-operative scrubsponge for use as a total body wash for patients prior to surgery

**Size**

Width of product	63 mm
Height of product	45 mm
Length of product	210 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-08

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**SURGISPONGE, impregnated with 4%  
chlorhexidine gluconate solution**

**Page 2**

**Material & colours**

Product material                      Polyurethane (PUR)  
Colour of product                      white

**Purity & certification**

Latex-free                                yes  
Pyrogen-free                            no  
Safety data sheet                        ja

**Packaging**

Minimum order qty.                    210  
Type of inner packaging                inner box  
Piece(s) / inner box                    35  
Piece(s) / case                         210  
Piece(s) / pallet                        3360  
Depth of case                            528 mm  
Width of case                            465 mm  
Height of case                           526 mm  
Case volume                             0,1283 cbm  
Weight of case                          8,2 kg  
EAN of inner box                        4038917150514  
EAN case                                 4038917036863

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-08

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**