



Product description

Order number	97.8720.500
Product description	SAHARA 4 100/240 V, for thawing and warming blood products in plastic bags before transfusion, max. capacity: 4 storage bags, operation via touchscreen, automatic aggregate status recognition by heating zone and ability to switch blood products conveniently while warming device is running

Product characteristics

Application	Warming
Power consumption during operation	1.000.00 W
Frequency	50 Hz - 60 Hz

Size

Width of product	574 mm
Height of product	348 mm
Length of product	554 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	Medical device
CE certified	CE 0482 (MedCert)
Batched	no

Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	1
Case volume	0.2095 cbm
Weight of product	23.7 kg
Weight of case	27.7 kg
EAN of inner box	4028339740127
EAN case	4028339740110

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

