



Product description

Order number 79.8710.575
Product description Paper roll, Consumables for protocol printer, white

Product characteristics

Application Warming

Size

Width of product 7.5 cm
Length of product 8.3 cm

Material & colours

Product material paper
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	10
Type of inner packaging	bag
Piece(s) / inner box	10
Piece(s) / case	50
Case volume	0.0319 cbm
Weight of product	0.2132 kg
Weight of case	11.86 kg
EAN of inner box	4028339753622
EAN case	4028339753615

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

