



Product description

Order number	97.8710.502
Product description	SAHARA-III 115 V, for thawing and warming blood products in plastic bags before transfusion, max. capacity: 3 storage bags, suitable for plasma, whole blood, erythrocyte concentrates and cryopreserved preparations

Product characteristics

Application	Warming
Supply voltage	115 V
Power consumption during operation	655 W
Frequency	50 Hz - 60 Hz

Size

Width of product	32 cm
Height of product	32,5 cm
Length of product	49,3 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	Medical device
CE certified	CE 0482 (MedCert)
Batched	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,129 cbm
Weight of product	14,64 kg
Weight of case	16,6 kg
EAN of inner box	4028339710229
EAN case	4028339710212

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

