



Product description

Order number	97.8710.550
Product description	Coding plug, for warming infusion solutions, tubes, instruments and contrast agent in the SAHARA-III

Product characteristics

Application	Warming
-------------	---------

Size

Width of product	4,5 cm
Length of product	6 cm

Purity & certification

Product category	Medical device
CE certified	CE
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0039 cbm
Weight of product	0,03 kg
Weight of case	0,22 kg
EAN of inner box	4028339753424
EAN case	4028339753417

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

