



Product description

Order number	79.8710.610
Product description	Replacement part, Gel compress for thawing out cryopreserved blood stem cell preparations in the SAHARA TSC

Product characteristics

Application	Warming
-------------	---------

Size

Width of product	16 cm
Length of product	23,5 cm

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE
Lot or serial number	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	5
Piece(s) / inner box	5
Piece(s) / case	5
Case volume	0,0041 cbm
Weight of product	0,45 kg
Weight of case	2,35 kg
EAN of inner box	4028339720426
EAN case	4028339720419

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

