



Product description

Order number	97.8710.600
Product description	Warming/thawing device, with active agitation plate, for thawing and warming blood products in plastic bags before transfusion, suitable for cryopreserved stem cells, incl. adaptation compress and warming tray, supply voltage: 230 V

Product characteristics

Application	Warming
Supply voltage	230 V
Power consumption during operation	655 W
Frequency	50 Hz - 60 Hz

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Width of product	32 cm
Height of product	32,5 cm
Length of product	49,3 cm

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE0482 (Medcert)
Lot or serial number	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,129 cbm
Weight of product	21,83 kg
Weight of case	23,79 kg
EAN of inner box	4028339720129
EAN case	4028339720112

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

