



Product description

Order number	97.8710.602
Product description	Warming/thawing device, with active agitation plate, for thawing and warming blood products in plastic bags before transfusion, suitable for cryopreserved stem cells, incl. adaptation compress and warming tray, supply voltage: 115 V

Product characteristics

Application	Warming
Supply voltage	115 V
Power consumption during operation	655 W
Frequency	50 Hz - 60 Hz

Size

Width of product	32 cm
Height of product	32,5 cm
Length of product	49,3 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	Medical device
CE certified	CE0482 (Medcert)
Batched	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,1639 cbm
Weight of product	21 kg
Weight of case	23,75 kg
EAN of inner box	4028339720228
EAN case	4028339720211

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

