



Product description

Order number	97.8710.620
Product description	Warming shell, Replacement part for storing stem cell preparations

Product characteristics

Application	Warming
-------------	---------

Size

Width of product	18 cm
Length of product	27 cm

Material & colours

Product material	aluminium
------------------	-----------

Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0096 cbm
Weight of product	2,45 kg
Weight of case	2,95 kg
EAN of inner box	4028339720525
EAN case	4028339720518

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

